



NON-RESORBABLE Ti-REINFORCED MEMBRANES

THE ORIGINAL, PATENTED, HIGH-DENSITY PTFE RE-INFORCED MEMBRANE (available in either 250µ or 150µ thickness)

- Membrane can be molded and shaped for tenting and space maintenance
- Regentex™ surface increases stability of membrane
- Dense PTFE prevents tissue ingrowth making removal of membrane easier when compared to removal of titanium mesh
- Impervious to bacteria (Data on file)
- Membranes can be trimmed to desired shape
- Titanium frame has little memory for easy placement



The patented Regentex™ surface helps stabilize the membrane, and the soft tissue flap.

CYTOPLAST™ Ti-Reinforce Anterior Narrow

	12mm x 24mm	
Ti250ANL-1	Ti150ANL-1	\$110
Ti250ANL-2	Ti150ANL-2	\$200



Ideal for narrow single-tooth extraction sites, especially where one or more bony walls are missing.

CYTOPLAST™ Ti-Reinforce Anterior Singles

	14mm x 24mm	
Ti250AS-1	Ti150ANL-1	\$140
Ti250AS-2	Ti150ANL-2	\$250



Ideal size and shape for single-tooth extraction sites, especially where one or more bony walls are missing.

CYTOPLAST™ Ti-Reinforce Ant/Post Perio

	13mm x 19mm	
Ti250PS-1	Ti150PS-1	\$180 Singles
Ti250PS-2	Ti150PS-2	\$325 Singles
Ti250PL-1	Ti150PL-1	\$200 Large
Ti250PL-2	Ti150PL-2	\$365 Large



Designed for grafting posterior extraction sites and limited ridge augmentation.

CYTOPLAST™ Ti-Reinforce Buccal

	17mm x 25mm	
Ti250BL-1	Ti150BL-1	\$180
Ti250BL-2	Ti150BL-2	\$325

Designed for grafting large buccal defects



CYTOPLAST™ Ti-Reinforce Posterior Singles

	20mm x 25mm	
Ti250PS-1	Ti150PS-1	\$180 Singles
Ti250PS-2	Ti150PS-2	\$325 Singles



CYTOPLAST™ Ti-Reinforce Posterior Large

	25mm x 30mm	
Ti250PS-1	Ti150PL-1	\$180 Singles
Ti250PS-2	Ti150PL-2	\$325 Singles



CYTOPLAST™ Ti-Reinforce Post Large T2

	25mm x 36mm	
Ti250PST2-1	Ti150PST2-1	\$215
Ti250PST2-2	Ti150PST2-2	\$390

Designed for grafting posterior extraction sites and limited ridge augmentation.



CYTOPLAST™ Ti-Reinforce XL / XLK

	30mm x 40mm	
Ti250XL-1	Ti150XL-1	\$250
Ti250XL-2	Ti150XL-2	\$450


Designed for grafting very large bony Defects especially ridge augmentation



CYTOPLAST™ Ti-Reinforce Ant Trans Crestal

NEW		38mm x 38mm	
Ti250ATC-1	Ti150ATC-1	\$230	
Ti250PTC-2	Ti150PTCL-2	\$420	

Designed for bony defects between teeth



CYTOPLAST™ Ti-Reinforce Posterior Distal

	38mm x 38mm	
Ti250PD-1	Ti150PD-1	\$280
Ti250PD-2	Ti150PD-2	\$500

CYTOPLAST™ Ti-Reinforce K2

	40mm x 50mm	
Ti250PD-1	Ti150PD-1	\$350
Ti250PD-2	Ti150PD-2	\$630

Designed for the largest bony defects, including augmentation



*** Products shown are NOT actual size.