

CYTOPLAST™ PTFE NEEDLE SIZING CHART

COMMON APPLICATIONS

- ⇒ Non-wicking monofilament construction
- ⇒ Soft and comfortable for patients
- ⇒ Excellent handling and knot security
- ⇒ Little or no package memory
- ⇒ Biologically inert



Benefits

- Non-resorbable medical grade PTFE
- Monofilament construction eliminates bacterial wicking
- Biologically inert & chemically non-reactive
- Not subject to weakening by tissue enzymes
- Excellent soft tissue response
- Soft; increased patient comfort
- Great handling, knot-tying & knot-holding characteristics



CS051819

New 300 series stainless steel needle - especially for dental implant and bone grafting procedures.

Size: USP 3-0

Needle: 19 mm 3/8 Circle Reverse Cutting

12 per box

\$100

CS0418

19 mm reverse cutting needle, which is popular for procedures where a longer needle is desired.

Size: USP 2-0

Needle: 19 mm 3/8 Circle Reverse Cutting

12 per box

\$100



CS0518

Most popular size especially for dental implant and bone grafting procedures.

Size: USP 3-0

Needle: 16 mm 3/8 Circle Reverse Cutting

12 per box

\$100

CS0618 RC

Most popular of the smaller CS-06 sutures, 16 mm reverse cutting needle. Often used for dental implant and bone grafting procedures.

Size: USP 4-0

Needle: 16 mm 3/8 Circle Reverse Cutting

12 per box

\$110

CS0618 PERIO

13 mm taper point needle. Often used on soft tissue grafts and suturing of delicate tissues that require an atraumatic needle.

Size: USP 4-0

Needle: 13 mm 1/2 Circle Taper Point

12 per box

\$110

CS0618 PREM

Offers over 50% longer cutting edges, giving needles a more slender, defined shape that greatly decreases the force applied to the needle during surgery.

Size: USP 4-0

Needle: 13 mm 3/8 Circle Premium Reverse Cutting

12 per box

\$110

*** Products shown are NOT actual size.



These sutures now available with black colored needles